



Design for a better life



TF Card



Air Filter **x5**



User Manual



Full Face Mask



Air Tubing



T-junction



Power Supply




Carrying Bag

# Auto BIPAP VM-7&VM-8



Patient-device synchrony  
Real-time displaying



Whisper quiet  
light and compact



Intelligent  
Breathing Algorithm




Efficient heating  
humidifier system

## Professional sleep study report

Unique sleep study report of whole week & month.

Easy operation, exquisite interface.



# Specification

Description	Auto BIPAP VM-7(ST25)	Auto BIPAP VM-8(ST30)
Work Mode	CPAP, Auto,S, T, ST, Auto S	CPAP, Auto,S, T, ST, Auto S, APCV
Sound	≤ 22 dBA @10cmH <sub>2</sub> O	
Screen	2.8 inch TFT, 320 x 240	
Pressure Range	4-25 cmH <sub>2</sub> O	4-30 cmH <sub>2</sub> O
IPAP	4-25 cmH <sub>2</sub> O	4-30 cmH <sub>2</sub> O
EPAP	4-20 cmH <sub>2</sub> O	4-20 cmH <sub>2</sub> O
Data Storage	> 5 years	
Power	Input: 100–240V~, 50–60Hz, 1.6A Max. Output: 24V = 2.5A	
Over-Temperature Protection	Stop automatically, alarm with 4 times beeps	
Parameter Lock	Yes	
Screen	Screen off in 15 minutes, press any button to highlight screen	
Backup Rate	5-50min, adjustable	
I/E	10-80%	
Ramp Time	0-60min, adjustable 1 min/Step	
Auto ON/OFF	Yes	
Inspirate Triggle	1-5 level	
Expirate Level	1-5 level	
Inspirate Sensitivity	1-5 level	
EPR(Expiratory Pressure Relief)	Yes, 0-3 level	
Auto-Altitude	Yes	
Leak Compensation	Yes	
Awake Detection	Yes	
Breath Event Detection	Yes	
Water Tank Capacity	200mL Max	
Humidifier Level	Integrated Heated, 0-5 level	
Alert	Power off, air leakage	
Dimensions(cm)	29L x 15W x 10H (machine)	
Carton Dimensions(cm)	34.5L x 20.5W x 25.5H (1pc carton)	
	39.5L x 34W x 25.5H (2pcs carton)	
	59.5L x 34W x 25.5H (3pcs carton)	
Gross Weight	4.5Kg (1pc/carton)	
	8.1Kg (2pcs/carton)	
	11.5Kg (3pcs/carton)	
Warranty	2 years	



**Hunan Ventmed Medical Technology Co., Ltd.**

Factory: Building No.9 , Area B of Shaoyang Electronic Information Industrial Park, No.18 Road, Shaoyang Economic Development Zone, Shaoyang City, Hunan, 422000, P.R. China

Tel.: +86-739-5188959 Web: www.ventmed.net